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General information about this medical device

The enclosed product is a custom-made product and was specifically designed, developed and manufactured according to the requirements of K05.002.FO Antrag für Sonderanfertigung / K05.004.FO Request for Custom Made Device DePuy Synthes - 41medical and K05.003.FO Bestellung Sonderanfertigung / K05.005.FO Order for Custom Made Device DePuy Synthes - 41medical. Accordingly 41medical cannot take the warranty for enclosed implant.

Possible side effects

- Implant failure due to wrong implant choice and / or overloading of the osteosynthesis.
- Allergic reactions due to material incompatibility
- Delayed healing due to vascular disorders
- Pain caused by the implant

Single-Use Products

Products intended for single use must not be re-used. Re-use or clinical processing (e.g. cleaning and re-sterilization) may compromise the structural integrity of the device and/or lead to device failure. This may result in patient injury, illness or death. Furthermore, re-use or clinical processing of single use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user. Do not reprocess soiled implants. Any 41medical implant that has been soiled by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though they may appear undamaged, the implants may have small defects and internal stress patterns that may cause material fatigue.

Non-Sterile Products

41medical products supplied in a non-sterile condition must be cleaned and steam sterilized prior to surgical use. Prior to cleaning, remove and dispose all original disposable packaging (e.g. silicone rubber guards, tip guards, protection caps, blisters, pouches, bags, packaging foam, card board etc.). Clean products before first and every use, and before returning for maintenance and repair. Prior to steam sterilization, place the product in an approved wrap or container. The first and most important step in decontaminating all re-usable instruments is thorough (manual and/or mechanical) cleaning and rinsing. Thorough cleaning is a complex process whose success depends on various interrelated factors: Water quality, quantity and type of cleaning agent, cleaning method (manual, ultrasonic bath, washer/disinfector), thorough rinsing and drying, proper product preparation, time, temperature, and thoroughness of the individual responsible for cleaning. Residual organic matter and/or a large number of microorganisms may reduce the effectiveness of the sterilization process.

Reprocessing of non-sterile 41medical implants and non-sterile single-use instruments

These recommendations are for processing of non-sterile 41medical implants and non-sterile single-use instruments. The information provided only for new and never before used implants and instruments. Explanted 41medical implants must not be reprocessed and disposed of according to the hospital protocol. Implants / instruments that have not yet been used but are contaminated must be disposed according to the hospital protocol. Contaminated implants must not be reprocessed.

Material Compatibility of 41medical implants and instruments in Clinical Processing

Stainless steels

– 41medical recommends disinfectants, cleaners or detergents with pH 7–11 for all stainless steels.

Titanium and its alloys

Titanium and titanium alloys are widely used as implant materials. On instruments titanium is used for only a few applications, mainly color coding of instruments. The surface of titanium alloys is also treated electrochemically (anodizing), producing a resistant oxide layer. Various color shades can be applied using this layer. Although anodized titanium and its alloys have good corrosion resistance, contact with strong alkaline detergents or disinfectants and solutions containing iodine or certain metal salts can lead to chemical attack and dissolution of the surface depending on the specific composition of the detergent. Consequently, 41medical recommends disinfectants, cleaners or detergents with pH 6 – 9.5. Products with a higher pH value, especially higher than pH 11, should only be used subject to the material compatibility requirements stated on the data sheet and other information from the manufacturer of the detergent.

Plastics

Various plastics are used for certain instrument parts, e.g. handles, radiolucent parts. All used plastics are able to withstand correct processing. Some of the plastics can become soft during steam sterilization, but do not undergo permanent deformation at normal sterilization temperatures below 140 °C. The material can, however, be damaged, for example by repeated immersion in disinfectants outside the pH range of 4 – 9.5 and by overstressing.

Recommended temperatures and pH levels

Material	Temperature	pH
Stainless steel	up to 149°C	7-11
Titanium alloys	up to 150°C	6-9.5
Plastics	up to 140°C	4-9.5

Important information

After use, all non-multiple-use 41medical instruments must be discarded according hospital protocol.

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Precautionary measures	<ul style="list-style-type: none"> After reception, check implants/instruments and verify identity, completeness, integrity and function. Before using any implants or instruments this cleaning and sterilization instruction has to be thoroughly read. DO NOT use steel-wool or any other detergent with abrasive properties. DO NOT use implants with a damaged surface. DO NOT use any cleaning solutions with iodine or high chlorine content. Products shall only be used in an ultrasound bath together with products that are made of metals and are of similar composition. ALL 41medical medical devices must be sterilized before use. The sterilization parameters are only valid for adequately pre-cleaned components. The parameters are only valid for properly installed, and adequately maintained and calibrated processing systems that comply with ISO 15883 and ISO 17665. Patients that are at risk for Creutzfeldt-Jacob-Disease (CJD) and its involved infections should be operated with single-use instruments. All instruments that have been used for a procedure in a patient with a suspected CJD exposure or confirmed diagnosis of CJD should be disposed after the procedure and the currently valid national regulations and recommendations have to be followed.
Limits of reprocessing	<ul style="list-style-type: none"> All custom made devices are single-use and must not be reprocessed.

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Processing instructions

First use at site of use	The non-sterile implants/instruments will be delivered in single product packaging, in pouches for medical use, clear packaging or similar products. On the packaging a product label is fitted. The products must not be used if the packaging is damaged. Sealings or protective films must be removed as well. Before use of any factory new implants they imperatively must have undergone a full processing cycle. Under no circumstances the cleaning process shall be omitted as during sterilization potential remnants or chemical residues on the instruments/implants (e.g. from packaging materials, maintenance fluids) could cause patches or coatings. The result of the cleaning process must be visually checked. The products must be deemed macroscopically clean.
Manual pre-cleaning	<ol style="list-style-type: none"> Clean the products with clear cold (< 40°C) water (drinking quality) until all visual contamination is removed. Adhering dirt should be removed with a soft brush (no wire or scratch brush). Cavities, cracks and slits should be flushed intensively for over 10 sec with a water-pistol or a similar device. Movable parts must be moved during the manual pre-cleaning. Clean the products in an ultrasound bath with an alkaline cleaning solution (0.5% neodisher® MediClean forte; solution/water temperature < 40°C). The sounding time is 10 min at a frequency of 35kHz. The instructions of the cleaning solution's manufacturer have to be followed. After the ultrasound-bath the products have to be re-cleaned with a brush. Difficult to access areas have to be flushed again with a water pistol or a similar device (> 10 sec). Rinse/flush the products for one minute with demineralized cold water. <p>Validated procedure: Equipment: Sink, soft brush, water pistol Cleaning solution: neodisher® MediClean forte Mixing ratio: 0.5% in normal tap water Temperature: Room temperature (23°C) Exposure time: 10 minutes</p>

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41medical AG