


# AODI Small External Fixator (AODI SEF)



THE AODI SEF is the perfect option to solve the setback of the instability of Kirschner wire fixation for fractures of the phalanx and the inapplicability of plates due to acute trauma and accessibility of soft tissues. In collaboration with the AO's Innovation Funding in Switzerland.

## Intended Use

The AODI Small External Fixator intended for temporary immobilization (a) in fractures of the metacarpal and phalanges of the hand for fracture stabilization and (b) in soft tissue surgery in the hand to reduce movement and thereby ensure undisturbed wound healing.

## Indications

- Fractures of the metacarpal and phalanges of the hand
- Soft tissue situations where temporary immobilization of the hand or finger is required

## Contraindications

- Advanced infection with incipient phlegmon
- Intolerance or hypersensitivity to the device materials
- Expected lack of patient compliance
- All types of fractures where adequate stability to bone healing by an external system is not guaranteed

Potential Side Effects and Risks please see instruction for use (IFU) on our website.

## Compatibility

No restrictions, no instruments needed for surgical procedure.

## Materials

- Carbon
- Polyurethane
- Stainless Steel

## Articles

- 100.126 AODI Small External Fixator 0.8x1.2mm
- 100.127 AODI Small External Fixator 1.0x1.2mm

In cooperation with:



# 41medical

Tannlimattweg 14 · 2544 Bettlach · Switzerland

T +41 32 645 41 41 · F +41 32 645 09 59 · info@41medical.com · www.41medical.com

Scan for more  
information  
41medical.com

