



Distal Radius Intramedullary Nail (DRIM-Nail)

The DRIM-Nail is a new, innovative treatment option for patients with extraarticular distal radius fractures, utilizing all benefits of a minimally invasive intramedullary fracture stabilization, resulting in a faster recovery, better and faster functional and cosmetic outcomes. In collaboration with DISRAD, Balgrist University Hospital and the AO's Innovation Funding in Switzerland.



Intended Use

The DRIM-Nail System is intended to be used for minimally invasive fixation of extra articular fractures of the distal radius.

Indications

Extra articular fractures of the distal radius (AO classification type 23-A2, 23-A3)

Contraindications

- Pre-existing or suspected infection at the fracture site
- Skeletally immature patients
- Insufficient bone quantity or quality
- Local or systemic use of substances that inhibit bone healing, such as but not only corticosteroids
- Any use outside of the defined indications

Potential Side Effects and Risks please see instruction for use (IFU) on our website.

Compatibility

- 41medical Locking Screw 2.5mm, L10-60mm
- 41medical DRIM-Nail Instrument System

Articles

- 100.130 DRIM-Nail
- 100.131 - 100.156 Locking Screw \varnothing 2.5mm L10-60mm

Materials

- Titanium Grade 2, ISO 5832-2, ASTM F67
- Stainless steel, ISO 5832-1, ASTM F138

In cooperation with:



41medical

Tannlimattweg 14 · 2544 Bettlach · Switzerland

T +41 32 645 41 41 · F +41 32 645 09 59 · info@41medical.com · www.41medical.com

Scan for more
information
41medical.com

