


Multi Axial Calcaneal Distractor (MACD)



MACD offers you essential support in surgical treatment of displaced calcaneal fractures. Distraction already helps significantly to reduce the fracture. The MACD offers additional axial rotation and Varus-Valgus correction, and stabilises your result for the follow-up steps of osteosynthesis. You can use it with your preferred surgical method (even ORIF). However, we would like to point out that it is developed and optimised for modern minimally invasive and mini-open procedures. In collaboration with PMU Innovations GmbH.

Intended Use

The MACD is intended to be used to anatomically reduce bones or bone fragments.

Indications

Calcaneal Fractures

Contraindications

- local or general prerequisites for surgery not given (local infection, ulcer, vascular occlusive disease)
- lack of compliance
- chronic alcohol and drug abuse

Potential Side Effects and Risks please see instruction for use (IFU) on our website.

Compatibility

Standard K-Wires with a diameter up to 3,2 mm

Articles

- 100.191 MACD Configuration I
- 100.192 MACD Configuration II

Materials

Stainless steel 1.4441 ASTM F138

In cooperation with:



41medical

Tannlimattweg 14 · 2544 Bettlach · Switzerland

T +41 32 645 41 41 · F +41 32 645 09 59 · info@41medical.com · www.41medical.com

Scan for more
information
41medical.com

